



NECTIN4 Antibody, HRP conjugated

Product Code	CSB-PA822274LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NY8
Immunogen	Recombinant Human Nectin-4 protein (91-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nectin-4 (Ig superfamily receptor LNIR) (Nectin cell adhesion molecule 4) (Poliovirus receptor-related protein 4) [Cleaved into: Processed poliovirus receptor-related protein 4], NECTIN4, LNIR PRR4 PVRL4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NECTIN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.