



WIPF2 Antibody, Biotin conjugated

Product Code	CSB-PA823478LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF74
Immunogen	Recombinant Human WAS/WASL-interacting protein family member 2 protein (18-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WAS/WASL-interacting protein family member 2 (WASP-interacting protein-related protein) (WIP- and CR16-homologous protein) (WIP-related protein), WIPF2, WICH WIRE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WIPF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.