



GATAD2B Antibody, Biotin conjugated

Product Code	CSB-PA823905LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI9
Immunogen	Recombinant Human Transcriptional repressor p66-beta protein (58-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional repressor p66-beta (GATA zinc finger domain-containing protein 2B) (p66/p68), GATAD2B, KIAA1150
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GATAD2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.