



RND1 Antibody, HRP conjugated

Product Code	CSB-PA835694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92730
Immunogen	Recombinant Human Rho-related GTP-binding protein Rho6 protein (1-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein Rho6 (Rho family GTPase 1) (Rnd1), RND1, RHO6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.