



SEMA4D Antibody, Biotin conjugated

Product Code	CSB-PA835707LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92854
Immunogen	Recombinant Human Semaphorin-4D protein (756-862AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-4D (A8) (BB18) (GR3) (CD antigen CD100), SEMA4D, C9orf164 CD100 SEMAJ
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SEMA4D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.