



SLC5A5 Antibody, Biotin conjugated

Product Code	CSB-PA835713LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92911
Immunogen	Recombinant Human Sodium/iodide cotransporter protein (547-643AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/iodide cotransporter (Na ⁺)/I ⁻ cotransporter) (Sodium-iodide symporter) (Na ⁺)/I ⁻ symporter) (Solute carrier family 5 member 5), SLC5A5, NIS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC5A5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.