



MOGAT1 Antibody, Biotin conjugated

Product Code	CSB-PA836275LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PD6
Immunogen	Recombinant Human 2-acylglycerol O-acyltransferase 1 protein (163-246AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-acylglycerol O-acyltransferase 1 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 1) (MGAT1) (Diacylglycerol O-acyltransferase candidate 2) (hDC2) (Diacylglycerol acyltransferase 2-like protein 1) (Monoacylglycerol O-acyltransferase 1), MOGAT1, DC2 DGAT2L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MOGAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.