



TSEN15 Antibody, FITC conjugated

Product Code	CSB-PA840989LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WW01
Immunogen	Recombinant Human tRNA-splicing endonuclease subunit Sen15 protein (1-171AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA-splicing endonuclease subunit Sen15 (SEN15 homolog) (HsSEN15) (tRNA-intron endonuclease Sen15), TSEN15, C1orf19 SEN15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TSEN15