



# SMG7 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA842128LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q92540
<b>Immunogen</b>	Recombinant Human Protein SMG7 protein (694-809AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein SMG7 (EST1-like protein C) (SMG-7 homolog) (hSMG-7), SMG7, C1orf16 EST1C KIAA0250
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	SMG7
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.