



SMG7 Antibody, Biotin conjugated

Product Code	CSB-PA842128LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92540
Immunogen	Recombinant Human Protein SMG7 protein (694-809AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein SMG7 (EST1-like protein C) (SMG-7 homolog) (hSMG-7), SMG7, C1orf16 EST1C KIAA0250
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SMG7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.