



SLC15A4 Antibody, Biotin conjugated

Product Code	CSB-PA843279LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N697
Immunogen	Recombinant Human Solute carrier family 15 member 4 protein (243-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 15 member 4 (Peptide transporter 4) (Peptide/histidine transporter 1) (hPHT1), SLC15A4, PHT1 PTR4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC15A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.