




SEMA6D Antibody

Product Code	CSB-PA847710LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFY4
Immunogen	Recombinant Human Semaphorin-6D protein (913-1068AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-6D, SEMA6D, KIAA1479
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SEMA6D

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, A549 whole cell lysate, 293 whole cell lysate

All lanes: SEMA6D antibody at 3µg/ml


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 120, 114, 112, 116, 118,

55, 68 kDa

Observed band size: 130 kDa



Immunofluorescence staining of HeLa cells with CSB-PA847710LA01HU at 1:133, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.