



WDFY2 Antibody, FITC conjugated

Product Code	CSB-PA850405LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P53
Immunogen	Recombinant Human WD repeat and FYVE domain-containing protein 2 protein (2-114AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat and FYVE domain-containing protein 2 (Propeller-FYVE protein) (Prof) (WD40- and FYVE domain-containing protein 2) (Zinc finger FYVE domain-containing protein 22), WDFY2, WDF2 ZFYVE22
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDFY2