



WDFY2 Antibody, Biotin conjugated

Product Code	CSB-PA850405LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P53
Immunogen	Recombinant Human WD repeat and FYVE domain-containing protein 2 protein (2-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat and FYVE domain-containing protein 2 (Propeller-FYVE protein) (Prof) (WD40- and FYVE domain-containing protein 2) (Zinc finger FYVE domain-containing protein 22), WDFY2, WDF2 ZFYVE22
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDFY2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.