



MRAP Antibody, HRP conjugated

Product Code	CSB-PA851545LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCY5
Immunogen	Recombinant Human Melanocortin-2 receptor accessory protein (59-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanocortin-2 receptor accessory protein (B27) (Fat cell-specific low molecular weight protein) (Fat tissue-specific low MW protein), MRAP, C21orf61 FALP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MRAP