



MRAP Antibody, HRP conjugated

Product Code	CSB-PA851545LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCY5
Immunogen	Recombinant Human Melanocortin-2 receptor accessory protein (59-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanocortin-2 receptor accessory protein (B27) (Fat cell-specific low molecular weight protein) (Fat tissue-specific low MW protein), MRAP, C21orf61 FALP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MRAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.