



ELP4 Antibody, FITC conjugated

Product Code	CSB-PA853408LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EB1
Immunogen	Recombinant Human Elongator complex protein 4 protein (117-223AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongator complex protein 4 (hELP4) (PAX6 neighbor gene protein), ELP4, C11orf19 PAXNEB
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ELP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.