



MCCC1 Antibody, Biotin conjugated

Product Code	CSB-PA853497ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RQ3
Immunogen	Recombinant Human Methylcrotonoyl-CoA carboxylase subunit alpha, mitochondrial protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylcrotonoyl-CoA carboxylase subunit alpha, mitochondrial (MCCase subunit alpha) (EC 6.4.1.4) (3-methylcrotonyl-CoA carboxylase 1) (3-methylcrotonyl-CoA carboxylase biotin-containing subunit) (3-methylcrotonyl-CoA:carbon dioxide ligase subunit alpha), MCCC1, MCCA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MCCC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.