



UBAC2 Antibody, HRP conjugated

Product Code	CSB-PA854093LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBM4
Immunogen	Recombinant Human Ubiquitin-associated domain-containing protein 2 protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-associated domain-containing protein 2 (UBA domain-containing protein 2) (Phosphoglycerate dehydrogenase-like protein 1), UBAC2, PHGDHL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBAC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.