



WNK3 Antibody, FITC conjugated

Product Code	CSB-PA861176LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYP7
Immunogen	Recombinant Human Serine/threonine-protein kinase WNK3 protein (1018-1181AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase WNK3 (EC 2.7.11.1) (Protein kinase lysine-deficient 3) (Protein kinase with no lysine 3), WNK3, KIAA1566 PRKWNK3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WNK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.