



SEMA6B Antibody, HRP conjugated

Product Code	CSB-PA862010LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3T3
Immunogen	Recombinant Human Semaphorin-6B protein (756-888AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-6B (Semaphorin-Z) (Sema Z), SEMA6B, SEMAN SEMAZ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SEMA6B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.