



# SLC19A3 Antibody, HRP conjugated

Product Code	CSB-PA866333LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZV2
Immunogen	Recombinant Human Thiamine transporter 2 protein (191-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiamine transporter 2 (ThTr-2) (ThTr2) (Solute carrier family 19 member 3), SLC19A3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC19A3