



PPA2 Antibody, HRP conjugated

Product Code	CSB-PA867137LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2U2
Immunogen	Recombinant Human Inorganic pyrophosphatase 2, mitochondrial protein (180-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inorganic pyrophosphatase 2, mitochondrial (EC 3.6.1.1) (Pyrophosphatase SID6-306) (Pyrophosphate phospho-hydrolase 2) (PPase 2), PPA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPA2