



HIPK2 Antibody, Biotin conjugated

Product Code	CSB-PA867138LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2X6
Immunogen	Recombinant Human Homeodomain-interacting protein kinase 2 protein (816-977AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-interacting protein kinase 2 (hHIPk2) (EC 2.7.11.1), HIPK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HIPK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.