



MBNL3 Antibody, HRP conjugated

Product Code	CSB-PA868322LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUK0
Immunogen	Recombinant Human Muscleblind-like protein 3 protein (116-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muscleblind-like protein 3 (Cys3His CCG1-required protein) (Muscleblind-like X-linked protein) (Protein HCHCR), MBNL3, CHCR MBLX39 MBXL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MBNL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.