



MBNL3 Antibody, FITC conjugated

Product Code	CSB-PA868322LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUK0
Immunogen	Recombinant Human Muscleblind-like protein 3 protein (116-211AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muscleblind-like protein 3 (Cys3His CCG1-required protein) (Muscleblind-like X-linked protein) (Protein HCHCR), MBNL3, CHCR MBLX39 MBXL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MBNL3