



MBNL3 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA868322LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NUK0 |
| Immunogen | Recombinant Human Muscleblind-like protein 3 protein (116-211AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Muscleblind-like protein 3 (Cys3His CCG1-required protein) (Muscleblind-like X-linked protein) (Protein HCHCR), MBNL3, CHCR MBLX39 MBXL |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | MBNL3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |