



INTS13 Antibody, Biotin conjugated

Product Code	CSB-PA868331LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVM9
Immunogen	Recombinant Human Integrator complex subunit 13 protein (573-706AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrator complex subunit 13 (Cell cycle regulator Mat89Bb homolog) (Germ cell tumor 1) (Protein asunder homolog) (Sarcoma antigen NY-SAR-95), INTS13, ASUN C12orf11 GCT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INTS13