



MOV10L1 Antibody, Biotin conjugated

Product Code	CSB-PA871604LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXT6
Immunogen	Recombinant Human RNA helicase Mov10l1 protein (336-425AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA helicase Mov10l1 (EC 3.6.4.13) (Moloney leukemia virus 10-like protein 1) (MOV10-like protein 1), MOV10L1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MOV10L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.