



# MEAF6 Antibody, Biotin conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA872530LD01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q9HAF1   |
| <b>Immunogen</b>           | Recombinant Human Chromatin modification-related protein MEAF6 protein (113-174AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Biotin   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4   |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Chromatin modification-related protein MEAF6 (MYST/Esa1-associated factor 6) (Esa1-associated factor 6 homolog) (Protein EAF6 homolog) (hEAF6) (Sarcoma antigen NY-SAR-91), MEAF6, C1orf149 CENP-28 EAF6 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling  |
| <b>Target Names</b>        | MEAF6  |