



# RXFP1 Antibody, Biotin conjugated

Product Code	CSB-PA875715LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBX9
Immunogen	Recombinant Human Relaxin receptor 1 protein (17-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Relaxin receptor 1 (Leucine-rich repeat-containing G-protein coupled receptor 7) (Relaxin family peptide receptor 1), RXFP1, LGR7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RXFP1