



LYRM4 Antibody, HRP conjugated

Product Code	CSB-PA875729LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD34
Immunogen	Recombinant Human LYR motif-containing protein 4 protein (1-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LYR motif-containing protein 4, LYRM4, C6orf149 ISD11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LYRM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.