



NDUFAF3 Antibody, Biotin conjugated

Product Code	CSB-PA880103LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BU61
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 3 protein (91-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 3, NDUFAF3, C3orf60
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NDUFAF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.