



SSBP3 Antibody, HRP conjugated

Product Code	CSB-PA880121LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWW4
Immunogen	Recombinant Human Single-stranded DNA-binding protein 3 protein (185-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-stranded DNA-binding protein 3 (Sequence-specific single-stranded-DNA-binding protein), SSBP3, SSDP SSDP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SSBP3