



SSBP3 Antibody, FITC conjugated

Product Code	CSB-PA880121LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWW4
Immunogen	Recombinant Human Single-stranded DNA-binding protein 3 protein (185-321AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-stranded DNA-binding protein 3 (Sequence-specific single-stranded-DNA-binding protein), SSBP3, SSDP SSDP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SSBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.