



MTMR12 Antibody, HRP conjugated

Product Code	CSB-PA880158LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0I1
Immunogen	Recombinant Human Myotubularin-related protein 12 protein (188-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 12 (Phosphatidylinositol 3 phosphate 3-phosphatase adapter subunit) (3-PAP) (3-phosphatase adapter protein), MTMR12, KIAA1682 PIP3AP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTMR12