



MTMR12 Antibody, Biotin conjugated

Product Code	CSB-PA880158LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C011
Immunogen	Recombinant Human Myotubularin-related protein 12 protein (188-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 12 (Phosphatidylinositol 3 phosphate 3-phosphatase adapter subunit) (3-PAP) (3-phosphatase adapter protein), MTMR12, KIAA1682 PIP3AP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTMR12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.