



SUN2 Antibody, Biotin conjugated

Product Code	CSB-PA883386LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UH99
Immunogen	Recombinant Human SUN domain-containing protein 2 protein (45-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SUN domain-containing protein 2 (Protein unc-84 homolog B) (Rab5-interacting protein) (Rab5IP) (Sad1/unc-84 protein-like 2), SUN2, FRIGG KIAA0668 RAB5IP UNC84B
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SUN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.