



HERC5 Antibody, HRP conjugated

Product Code	CSB-PA883435LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UII4
Immunogen	Recombinant Human E3 ISG15--protein ligase HERC5 protein (153-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ISG15--protein ligase HERC5 (EC 2.3.2.-) (Cyclin-E-binding protein 1) (HECT domain and RCC1-like domain-containing protein 5), HERC5, CEB1 CEBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HERC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.