



PIGQ Antibody, FITC conjugated

Product Code	CSB-PA883583LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRB3
Immunogen	Recombinant Human Phosphatidylinositol N-acetylglucosaminyltransferase subunit Q protein (616-760AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol N-acetylglucosaminyltransferase subunit Q (EC 2.4.1.198) (N-acetylglucosamyl transferase component GPI1) (Phosphatidylinositol-glycan biosynthesis class Q protein) (PIG-Q), PIGQ, GPI1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PIGQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.