



SDHAF2 Antibody, Biotin conjugated

Product Code	CSB-PA885749LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX18
Immunogen	Recombinant Human Succinate dehydrogenase assembly factor 2, mitochondrial protein (30-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Succinate dehydrogenase assembly factor 2, mitochondrial (SDH assembly factor 2) (SDHAF2), SDHAF2, C11orf79 PGL2 SDH5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SDHAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.