



NFKBIL1 Antibody, Biotin conjugated

Product Code	CSB-PA887001LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBC1
Immunogen	Recombinant Human NF-kappa-B inhibitor-like protein 1 protein (185-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NF-kappa-B inhibitor-like protein 1 (Inhibitor of kappa B-like protein) (I-kappa-B-like protein) (IkappaBL) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 1), NFKBIL1, IKBL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NFKBIL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.