



IFT172 Antibody, HRP conjugated

Product Code	CSB-PA887019LB01HU
Abbreviation	Intraflagellar transport protein 172 homolog
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UG01
Immunogen	Recombinant Human Intraflagellar transport protein 172 homolog protein (1368-1502AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Required for the maintenance and formation of cilia. Plays an indirect role in hedgehog (Hh) signaling, cilia being required for all activity of the hedgehog pathway (By similarity).
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 172 homolog, IFT172, KIAA1179
Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IFT172