



IFT172 Antibody, Biotin conjugated

Product Code	CSB-PA887019LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UG01
Immunogen	Recombinant Human Intraflagellar transport protein 172 homolog protein (1368-1502AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 172 homolog, IFT172, KIAA1179
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IFT172
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.