



GPR174 Antibody, Biotin conjugated

Product Code	CSB-PA887146LD01HU
Abbreviation	Probable G-protein coupled receptor 174
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXC1
Immunogen	Recombinant Human Probable G-protein coupled receptor 174 protein (291-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Putative receptor for purines coupled to G-proteins.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 174, GPR174
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR174