



HES7 Antibody, Biotin conjugated

Product Code	CSB-PA887160LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYE0
Immunogen	Recombinant Human Transcription factor HES-7 protein (127-192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor HES-7 (hHes7) (Class B basic helix-loop-helix protein 37) (bHLHb37) (Hairy and enhancer of split 7) (bHLH factor Hes7), HES7, BHLHB37
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HES7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.