



REG4 Antibody, HRP conjugated

Product Code	CSB-PA887169LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYZ8
Immunogen	Recombinant Human Regenerating islet-derived protein 4 protein (27-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regenerating islet-derived protein 4 (REG-4) (Gastrointestinal secretory protein) (REG-like protein) (Regenerating islet-derived protein IV) (Reg IV), REG4, GISP RELP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	REG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.