



REG4 Antibody, HRP conjugated

Product Code	CSB-PA887169LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYZ8
Immunogen	Recombinant Human Regenerating islet-derived protein 4 protein (27-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regenerating islet-derived protein 4 (REG-4) (Gastrointestinal secretory protein) (REG-like protein) (Regenerating islet-derived protein IV) (Reg IV), REG4, GISP RELP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	REG4