



FAM234A Antibody, Biotin conjugated

Product Code	CSB-PA887949LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0X4
Immunogen	Recombinant Human Protein FAM234A protein (115-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM234A (Protein ITFG3), FAM234A, C16orf9 ITFG3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM234A