



FAM234A Antibody, Biotin conjugated

Product Code	CSB-PA887949LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0X4
Immunogen	Recombinant Human Protein FAM234A protein (115-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM234A (Protein ITFG3), FAM234A, C16orf9 ITFG3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM234A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.